



Electrostimulator muscular MD6150

- electrostimuloarele musculare LAICA functioneaza printr-o stimulare electrica de frecventa joasa care reproduce mecanismul natural al contractiilor musculare
- usor de folosit si sigur, non-toxic, aparatul pentru electrostimulare nu strapunge pielea
- este ideal a se folosi acasa
- construit cu ecran mare LCD usor de citit
- acest electrostimulator LAICA ofera electrostimulare musculara si electrostimulare nervoasa transcutanata
- 2 canale
- 10 programe
- 4 electrozi
- 15 nivele de intensitate
- 3 baterii alcaline (tip 1.5V AAA) incluse
- design italian
- produs construit in conformitate cu directiva 93/42/EEC pentru "dispozitive medicale"
- 2 ani garantie

Detalii tehnice:

Dimensiuni: L=1,5 Lat=7 H=16,7 (cm)

Greutate: 0,1 (kg)


Dimensiunile ecranului LCD: 2 x 4 (cm)

Codul de bare: 8013240401454

Ambalaj final: 4 produse


Dimensiunile ambalajului: L=25 Lat=35,5 H=36,5 (cm)

Greutatea ambalajului: 3,3 (kg)



Culoare: negru / argintiu



CE 0434 (dir. 93/42/EEC) Italian Design 

SC ECOMEDIS SRL
Sos. Borsului, nr 35, ap 15, et 3
Oradea, jud. Bihor
Tel. +40 0259464748
Fax. +40 0259464749
E-mail: office@laica.ro
Web site: www.laica.ro